

READ THESE INSTRUCTIONS

Coren[®]PS PRACTICE CARD



USING EARPIERCER
with REUSABLE
FLEXHANDLE



*The COREN[®]PS Earpiercer
can be used by itself or
loaded into the
Reusable FLEXHANDLE.*

See reverse side for instructions
on using the COREN[®]PS without
the Flexhandle

SIMPLE PROCEDURE...

- Grip FLEXHANDLE securely as shown above
- Close gently at spot of piercing
- Slowly increase pressure to break seal - instantly piercing lobe
- Release pressure (open handle)
- Slide COREN[®]PS downward to disengage it from clasp (See Fig. 4 on reverse side)

For frequently asked questions, or ordering information,
please visit www.CorenPS.com or call 1-877-267-3677.

THE EARPIERCER SHAPE IS A TRADEMARK OF HASTINGS CORP • U.S. AND FOREIGN PATENTS



Practice
HERE

INSTRUCTIONS FOR USING COREN®PS EARPIERCER *without* FLEXHANDLE




Figure 1



Figure 2



Figure 3




Figure 4

SIMPLE PROCEDURE...

- Remove the COREN®PS from its sterile package
- Hold securely as shown (Fig. 1)
- Close gently at spot of piercing (Fig. 2)
- Slowly increase pressure to break seal - instantly piercing lobe (Fig. 3)
- Release pressure
- Slide COREN®PS downward to disengage it from clasp (Fig. 4)

See reverse side for practicing, and instructions on using the COREN®PS with the Flexhandle